

Process:
Subprocess:
Product Family:

FOAMING
API
MEGATRON®



PHARMACEUTICAL

Foamed patches for efficient substance release



THE APPLICATION SCOPE / PURPOSE

Nicotine patches are the most commonly used items to fight the addiction of smoking. The release of nicotine should be as efficient as possible. The surface in contact with the skin can be enlarged by tiny drops of emulsions or foam bubbles in the patch generated through dispersion solutions from Kinematica. This guarantees an acceleration of the nicotine release reaction.

Get full access to the Application Note by clicking this link:

GET FULL VERSION

Kinematica AG

Werkstrasse 7 c-d / CH-6102 Malters / Switzerland
Phone +41 41 259 65 65 / Fax +41 41 259 65 75
info@kinematica.ch / kinematica.ch

🇨🇭 Swiss made dispersing solutions for you. Since 1962.